## Laboratory Sample Handling



## Pipette Aid Safe, Easy and Economical Pipetting

## LAB-378860000

This pipette aid is used by squeezing, then releasing the bulb to draw fluid into the pipette. The pipette aid is then removed, and the pipette top quickly covered with a gloved finger to control and release fluid. The bulb is Vikem® vinyl with a tapered polyethylene chuck for quick removal when pipette is filled. 6 per pack.

**Dimensions:** 5 dia x 10 cm H (1.9 x 3.9")



| To Order      |             |
|---------------|-------------|
| Model No.     | Description |
| LAB-378860000 | Pipette aid |